


How the Acne LED Light Works

LED blue light therapy helps to kill the acne causing bacteria, reduce inflammation and provide a general rejuvenation effect on the skin. Usually, in the case of inflamed acne, topical products will bring down the inflammation rather quickly. However, there are times when the skin needs a little more help in bringing down inflammation. This is where the LED can help.



The combination of blue, red and/or infrared lights in the LightStim Acne Light get down to the core of the acne causing bacteria. The red lights help to reduce inflammation while at the same time, the blue lights help to kill the acne causing bacteria.

HOW YOU WILL USE IT

1. To cleanse the top of the LED, wipedown the lights down with a cotton round soaked or sprayed with alcohol.
2. Cleanse your face and dry it
3. Plug in the LED and turn on the switch
4. Place the LED right on the skin in the area that has inflammation
5. Leave it on each area for 3 minutes
6. When you have covered the whole area of inflammation, then do the rest of your skincare routine

This is not a substitute for products; nor is it a substitute for your scheduled appointments with us.

****You still need to use your products, come in for regular treatments and get your routine adjusted every 2 weeks. This device must be returned to us within 60 days. If the LED is not returned, your credit card on file will be charged.**

I agree to return the LightStim LED by _____

This content is being provided for informational purposes only. This is not meant to be a substitute for professional medical advice, diagnosis or treatment. Face Reality is not offering, recommending, or endorsing any specific dietary guidelines, prescription drug, pharmacy or other information on the site. Face Reality provides no warranty for any information. Please seek medical advice before starting, changing or terminating any medical treatment.